ENGLISH
INSTRUCTIONS FOR USE • I.B.S™ Osteosynthesis screws STERILE IMPLANTS FOR OSTEOSYNTHESIS OR ORTHOPEDIC SURGERY SINGLE USE

In accordance with the directive 93/42/EEC relative to medical devices and its amendments, this product must be handled and/or implanted by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.

1. Description of the medical device

The I.B.S.™ screws are:

osteosynthesis bone screws, existing in different models, diameters and lengths
The Implants are manufactured in:
- titanium alloy in accordance with the standard NF ISO 5832-3 and ASTM F136
The medical devices are sold sterile.
Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation.
These devices do not contain phtalates unless this is indicated on the label.

2. Indications

The I.B.S.™ compression and neutralization osteosynthesis screws are intended for:
- The fixation of arthrodiesis, osteotomies or fractures of long or short bones of the upper and lower limbs
- Osteosynthesis requiring a mono or bicortical compression

The size of the chosen screw should be adapted to the specific indications.

The I.B.S.™ snap off screws are intended for:
- Bone reconstruction, osteotomy and fracture fixation of bones appropriate for the size of the device
Examples include:
Mono-cortical fixation of small bone fragments
Weil osteotomy
Osteotomies and fractures fixation in the foot and hand.

3. Contraindications

The implant should not be used in a patient who has currently, or who has history of:

- systemic, acute or chronic inflammations, .
- active infections,
- sensitivity/allergies to the implant materials (cf paragraph 1).

4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic degenerative, acute osteopenia, bone tumor, etc...)
- Use of steroid derivatives, chemotherapy,... etc.
- Drug and/or alcohol and/or smoke addiction and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions

Complications may include but are not limited to:

- Pain, discomfort or abnormal sensations due to presence of the implant

- Bending, loosening and/or breakage
- Migration of the implant position
- Bone loss due to stress shielding
- Deformation recurrence, loss of correction
- Delayed union or pseudarthrosis
- Infections, hematoma, allergy, thrombosis

5. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon.

Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

5.1. Preoperative

- Read carefully the Surgical Technique.
- Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.
- Prepare all implants and instruments necessary for the surgical procedure. Do not attempt a surgical procedure with non-functional, broken, suspect or damaged instrument.

- The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.

- The size and number of implants should be assessed based on the preoperative X-rays.
- After measurements, some implants from each adequate size should be made available to have a sufficient inventory for the surgery.

- An additional implant from each size should be made available to replace any implant that might be accidentally contaminated during the surgery.

- Before the first surgery, the surgeon and assistants should manipulate the instruments to familiarize themselves with the material.

5.2. Peroperative

The surgery should be performed by a surgeon with adequate background in orthopaedics and with respect to the different steps described in the Surgical Technique.

- Implants should be handled with care to avoid any scratch (risk of incipient break).
- Under no circumstances should the implant be modified.
- Alternate fixation methods should be available intraoperatively.
- Opening of the instruments set must be done according to aseptc condition.
- Check packaging and labeling integrity before use.
- Do not use any implant for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments).

5.3. Postoperative

- The patient should be advised that a second more minor procedure for the removal of the implants may be necessary.
- Some X-rays should be periodically done to check the postoperative progress and prevent any complication.

- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing or excessive muscular activity. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until stability is evident.

- The patient should be encouraged to report to his/her surgeon any unusual changes of the operated extremity. If evidence suggests loosening of the implant (particular pain and progressive changes in the radiographs) an intensified schedule of check-ups is advised and new warning and instructions to the patient may be appropriate regarding further activity restrictions.

- The patient should be encouraged to receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.

5.4. Re-use / Re-sterilization

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use device may create risks of contamination from one patient to another or the use. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even though they may appear undamaged, the implants may have defects and internal stress patterns that may cause material fatigue.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

6. Removal of the implant after healing

Particularly in young active patients, implants may loosen, fracture, migrate, increase the risk of infection, cause pain, or stress shield bone – even after normal healing. The surgeon should consider the risks and benefits when deciding whether or not to remove an implant. Removal should be followed by careful postoperative management to avoid re-fractures. If the patient is older and has a low activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of another surgery.

7. Storage

Store in dry place.

8. MRI/SCANNER

The patient should be asked to systematically mention that he/she was implanted with a surgical orthopedic / osteosynthesis implant.

The medical device(s) has not been evaluated for safety and compatibility in the MR environment. The medical device(s) has not been tested for heating or migration in the MR environment.

9. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

WARNING (USA): This device has not been approved for fixation or maintaining purposes in the back part of the intervertebral body (neck of vertebra) or at cervical, thoracic or lumber levels.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

INFORMATION

Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer.

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